

ORDER FORM

CORZINO Basic

Release date 2025/10



Company	Invoice address	Order number / reference	
Street		Customer reference	
E-Mail		Date	
City, Zip code		Offer	
Telephone		Order	
Company	Delivery address		
Street	City, Zip code		

Standard equipment:	Size / order codes			
• 2-piece footrest height adjustable and foldable • backrest angle adjustable (2 positions) • 5-point immmo. harness with crotch belt and shoulder padding • backrest height extension (size 38, 42 only) • seat tension adjustment • 175/265 mm solid wheels	30	34	38	42
	front swivel castor-wheels /FSCW			
	CNB30RP02	CNB34RP02	CNB38RP02	CNB42RP02
	€	€	€	€



Seat width [mm]	300	340	380	420
Seat depth [mm]	280	330	380	380
Backrest height [mm]	550	650	650 (820)	650 (830)
Seat to footplate [mm]	290-420	300-440	300-440	330-490
Backrest angle	91°,93°	91°,94°	90°,93°	90°,93°
Seat angle	22°,21°	23°,22°	22°,21°	23°,22°
Maximum load	45,0 kg	55,0 kg	75,0 kg	75,0 kg

colour frame / pattern colour	
frame	K1 - Anthracite
fabric patterns	G02 Tekno ruby red G03 Tekno dark blue G06 Tekno grass green
	G07 Tekno black

CORZINO Accessories



Section E: Positioning and fixation - Head

No.	Desc. / Codes	CNB 30	CNB 34	CNB 38	CNB 42	
E1	Headrest CLASSIC RPBE001	UNI	UNI	UNI	UNI	€
E2	Headwings BASIC RPBE002	M	M	L	L	€





Section G: Positioning and fixation - Trunk

G1	2-Point Safety Belt RPBG001	S	M	M	M	€
G2	SP - 5-Point Immo. Belt CLASSIC RPBG003	S(SE)	M(SE)	M(SE)	M(SE)	€
G3	5-Point immo.belt COMFORT RPBG006	S	M	M	M	€
G4	H-Style Harness RPBG008	S	M	M	L	€
G5	Immobilisation Waistcoat RPBG009	UNI	UNI	UNI	UNI	€
G6	Chest immo. belt RPBG010	S	M	M	L	€
G7	Extended side protection RPBG011	S	M	L	L	€






Section I: Positionig and fixation - Pelvis / Abduction

No.	Desc. / Codes	CNB 30	CNB 34	CNB 38	CNB 42	
I1	Comfort side protection RPBI001	S	M	M	M	€
I2	Abduction block CLASSIC RPBI003	S	M	M	M	€







Section K: Positionig and fixation - Arms

K1	Armrest foldable, (set) RPBK001	M	M	L	L	€
K2	Work board (soft), CLASSIC RPBK101	UNI	UNI	UNI	UNI	€



Section M: Positionig and fixation - Legs



M1	Thigh fixation, CLASSIC RPBM002	XS/S	S	M	M	€
----	------------------------------------	------	---	---	---	---



Section P: Protection and comfort





P1	Canopy buggy BASIC RPBP003	S	M	L	L	€
P2	Rain cover RPBP051	M	M	L	L	€
P3	Footmuff - summer CLASSIC RPBP101	M	M	L	L	€

P4	Footmuff - winter RPBP151	S	M	M	L	€
P5	Seat cover COMFORT RPBP201	M	M	L	L	€







Section R: Bags, baskets

R1	Basket - textile, RPBR001	M	M	L	L	€
----	------------------------------	---	---	---	---	---




Section T: Accessory - specialities

T1	Lead antist.(ground.)- RPBT005	UNI	UNI	UNI	UNI	€
----	-----------------------------------	-----	-----	-----	-----	---






Section U: Accessory - sets

U1	Footrest 3D-Flexi, 2-piece (2PFF) (size:30) RPBU00130	UNI (SE)			€	
U2	Footrest 3D-Flexi, 2-piece (2PFF) (size:34,38,42) RPBU00134		UNI (SE)	UNI (SE)	UNI (SE)	€
U3	Drum-brake adapt.set (w. hand control), (38,42) RPBU10138			UNI	UNI	€



Section V: Transportation / anchoring

V1	'TIE DOWN' set (EU markets only) RPBV001	UNI	UNI	UNI	UNI	€
V2	'TIE DOWN' set (USA, CDN markets only) RPBV002	UNI	UNI	UNI	UNI	€



Sizes:

XS - extra small
 S - small
 M - medium
 L - large
 XL - extra large
 UNI - universal size (one size only)

Other notes:

(SE) Standard Equipment
 (empty position) - accessory is not available