


# ORDER FORM

## CORZINO Basic













Release date 2024/11



Company	Invoice address	Order number / reference	
Street		Customer reference	
E-Mail		Date	
City, Zip code		Offer	
Telephone		Order	
Company	Delivery address		
Street	City, Zip code		

Standard equipment:	Size / order codes			
	30	34	38	42
<ul style="list-style-type: none"> <li>• 2-piece footrest height adjustable and foldable</li> <li>• backrest angle adjustable (2 positions)</li> <li>• 5-point immmo. harness with crotch belt and shoulder padding</li> <li>• backrest height extension (size 38, 42 only)</li> <li>• seat tension adjustment</li> <li>• 175/265 mm solid wheels</li> </ul>	front swivel castor-wheels /FSCW			
	CNB30RP02	CNB34RP02	CNB38RP02	CNB42RP02
	€	€	€	€

Seat width [mm]	300	340	380	420
Seat depth [mm]	280	330	380	380
Backrest height [mm]	550	650	650 (820)	650 (830)
Seat to footplate [mm]	290-420	300-440	300-440	330-490
Backrest angle	91°,93°	91°,94°	90°,93°	90°,93°
Seat angle	22°,21°	23°,22°	22°,21°	23°,22°
Maximum load	45,0 kg	55,0 kg	75,0 kg	75,0 kg



colour frame / pattern colour							
<b>frame</b>	<b>K1 - Anthracite</b>						
<b>fabric patterns</b>	<table style="width: 100%; text-align: center;"> <tr> <td>G02  Tekno ruby red</td> <td>G03  Tekno dark blue</td> <td>G06  Tekno grass green</td> </tr> <tr> <td>G07  Tekno black</td> <td></td> <td></td> </tr> </table>	G02  Tekno ruby red	G03  Tekno dark blue	G06  Tekno grass green	G07  Tekno black		
	G02  Tekno ruby red	G03  Tekno dark blue	G06  Tekno grass green				
G07  Tekno black							

## CORZINO Basic Accessories

### Positioning and Fixation




#### Head

No.	Desc. / Codes	CNB 30	CNB 34	CNB 38	CNB 42	
A1	Headrest adjustable RPRB002	UNI	UNI	UNI	UNI	€
A2	Headwings BASIC RPRB013	M	M	L	L	€
A3	Backrest Height Extension			UNI(SE)	UNI(SE)	€




#### Trunk

B1	SP - 5-Point Immo. Belt CLASSIC RPRB010	S(SE)	M(SE)	M(SE)	M(SE)	€
B2	Chest immo. belt RPRB004	S	M	M	L	€
B3	Immobilisation Waistcoat RPRB005	UNI	UNI	UNI	UNI	€
B4	H-Style Harness RPRB007	S	M	M	M	€
B5	Comfort side protection RPRB020	S	M	M	M	€
B6	Extended side protection RPRB032	S	M	L	L	€









## Pelvis / abduction

No.	Desc. / Codes	CNB 30	CNB 34	CNB 38	CNB 42	
C1	Abduction Block RPRB015	S	M	M	M	€
C2	2-Point Safety Belt RPRB019	S	M	M	M	€
C3	Thigh fixation RPRB024	S	S	M	M	€






## Legs / Arms

D1	Footrest BASIC, 2-piece; RPRB04130	UNI(SE)				€
D2	Footrest BASIC, 2-piece RPRB04134		UNI(SE)	UNI(SE)	UNI(SE)	€
D3	Work board (soft) RPRB018	UNI	UNI	UNI	UNI	€
D4	Armrest foldable (set) RPRB02302	UNI	UNI			€
D5	Armrest foldable (set) RPRB02304			UNI	UNI	€



## Protection / transport accessories

No.	Desc. / Codes	CNB 30	CNB 34	CNB 38	CNB 42	
E1	Canopy buggy BASIC RPRB00102	UNI	UNI	UNI	UNI	€
E2	Rain cover - M RPRB01701	M	M			€
E3	Rain cover - L RPRB01702			L	L	€
E4	COMFORT Seat Cover RPRB003	M	M	L	L	€
E5	Basket - textile RPRB00801	UNI	UNI			€
E6	Basket - textile RPRB00802			UNI	UNI	€
E7	Footmuff - Summer RPRB012	M	M	L	L	€
E8	Footmuff - Winter RPRB025	S	M	M	L	€



### Sizes:

XS - extra small  
 S - small  
 M - medium  
 L - large  
 XL - extra large  
 UNI - universal size (one size only)

### Other notes:

(SE) Standard Equipment  
 (empty position) - accessory is not available