

ORDER FORM

PIPER Comfort

Release date 2025/10



Company	Invoice address	Order number / reference	
Street		Customer reference	
E-Mail		Date	
City, Zip code		Offer	
Telephone		Order	
Company	Delivery address		
Street	City, Zip code		

Standard equipment:	Size / order codes
<ul style="list-style-type: none"> • seat depth non-adjustable • extended backrest height non-adjustable • backrest recline adjustable • footrest height adjustable, detachable • seat width adjustable by FLEXI-pads set (optional eq.) • 5-point immo.harness with crotch belt and shoulder padding • 175/225 mm solid wheels, front with swivel lock • detachable wheels, forks • one-piece foldable handle • basket, transportation lock 	front swivel castor-wheels /FSCW
	<p>PIF34RP01</p> <p>€</p>




Seat width [mm]	260/step/340	with Seat inlay / footrest raiser	230/270
Seat depth [mm]	290		240
Backrest height [mm]	680		680
Seat to footplate [mm]	240-370		200/230
Backrest angle	90° - 125°		90° - 125°
Seat angle	25°		25°
Maximum load	50,0 kg		50,0 kg

colour frame / pattern colour	
frame	K2 - Anthracite
fabric patterns	<p>G17  Denim blue</p> <p>G20  Teal green</p> <p>G21  Grey-black</p>
	<p>G22  Antique pink</p> <p>G25  Alphabet grey</p>

PIPER Comfort Accessories


Section E: Positionig and fixation - Head

No.	Desc. / Codes	PIF34	
E1	Headrest buggy adj. CLASSIC RPBE003	UNI	€





Section G: Positionig and fixation - Trunk

G1	SP-Immo.belt 5-point CLASSIC RPBG004	UNI(SE)	€
G2	H-style harness CLASSIC RPBG008	(S (1)) M	€
G3	Seat inlay (seat size reduction) RPBG051	UNI	€




Section I: Positionig and fixation - Pelvis / Abduction

I1	FLEXI-pad set (adj. side cushion set) RPBI002	UNI	€
I2	Abduction block CLASSIC RPBI003	(S (1)) M	€




Section M: Positionig and fixation - Legs

M1	Feet fixation- pair RPBM001	UNI	€
M2	Thigh fixation, CLASSIC RPBM002	(XS (1)) S	€
M3	Footplate booster ((soft, 5+5cm height) RPBM051	UNI	€






Section P: Protection and comfort

P1	Canopy buggy CLASSIC RPRB02901	UNI	€
P2	Rain cover buggy RPBP051	M	€
P3	Footmuff PIPER - Summer (extend.) RPBP102	UNI	€



Section V: Transportation / anchoring

V1	TIE DOWN' set PIF (EU markets only) RPBV003	UNI(1)	€
V2	TIE DOWN' set PIF (USA,CDN markets only) RPBV004	UNI(2)	€



Sizes:

S - small
M - medium
UNI - universal size (one size only)

Other notes:

(SE) Standard Equipment
(1) recomended with seat inlay or for smaller children