

ORDER FORM

TOM4 Classic





Release date 2018/05



| | | | |
|----------------|------------------|--------------------------|--|
| Company | Invoice address | Order number / reference | |
| Street | | Customer reference | |
| E-Mail | | Date | |
| City, Zip code | | Offer | |
| Telephone | | Order | |
| Company | Delivery address | | |
| Street | City, Zip code | | |

| Standard equipment: | Size / order codes | | | |
|---|----------------------------------|-----------|----------|-----------|
| | MINI | STD | MAXI | SM 42 |
| <ul style="list-style-type: none"> • detachable and reversible seat-unit • height adjustable backrest and depth adjustment of the seat-base. • width adjustable seat by FlexiPads set (optional equipment) • seat with fully adjustable tilt in space function • fully reclining backrest • angle and height adjustable footrest • height adjustable ergonomic handlebar • 200/300 mm solid wheels • detachable soft side-cushions • detachable grab bar with cover | front swivel castor-wheels /FSCW | | | |
| | T4CWKPMYY | T4CWKPYYY | T4CWKPXY | T4CWKPS2Y |
| | € | € | € | € |

| | | | | |
|-----------------------|-----------|-------------|-------------|-----------|
| Backrest height [mm] | 450-600 | 520-680 | 570-680 | 570-680 |
| Seat depth [mm] | 270-310 | 270-310 | 320-430 | 320-430 |
| Seat width [mm] | 260/300 | 240/280/340 | 280/320/380 | 425 |
| Seat to footrest [mm] | 140-280 | 200-350 | 240-400 | 240-400 |
| Backrest angle | 0° - 175° | 0° - 175° | 0° - 165° | 0° - 165° |
| Seat angle | +/-27° | +/-27° | +/-25° | +/-25° |
| Maximum load | 35,0 kg | 35,0 kg | 50,0 kg | 75,0 kg |

| colour frame / pattern colour | |
|-------------------------------|--|
| frame | K2 - Anthracite |
| fabric patterns | <p>R43  Dark lime green</p> <p>R44  Lavender blue</p> <p>R45  Ruby red</p> <p>R50  Pea green</p> |

TOM4 Accessories

Positioning and Fixation

Head



| No. | Desc. / Codes | T4C MINI | T4C STD | T4C MAXI | T4C SM42 | |
|-----|------------------------------|----------|---------|----------|----------|---|
| A1 | Headrest CLASSIC RPRK002 | S | M | M | M | € |
| A2 | Headwings CLASSIC RPRK049 | S | M | M | M | € |
| A3 | Neckrest CLASSIC RPRK028 | S | M | M | M | € |



Trunk


| | | | | | | |
|-----|---|---------|-----|-----|-----|---|
| B1 | Side trunkrest CLASSIC RPRK003 | S | M | M | L | € |
| B2 | Chest immo. belt CLASSIC RPRK012 | S | M | M | L | € |
| B3 | Hip rest CLASSIC RPRK006 | M | M | L | L | € |
| B4 | Lumbar support W25 RPRK04302 | UNI | UNI | | | € |
| B5 | Lumbar support W35 RPRK04303 | | | UNI | UNI | € |
| B6 | Backrest wedge (reduction seat depth 245 mm) RPRK081 | S | | | | € |
| B7 | 5-point immobilisation belt CLASSIC RPRK047 | S | M | M | L | € |
| B8 | 4-point immobilisation belt CLASSIC RPRK004 | S | M | M | L | € |
| B9 | Chest-belt addition to 4P-belt RPRK041UN | UNI | UNI | UNI | UNI | € |
| B10 | Immo. Waistcoat CLASSIC RPRK009 | (XS) S | M | M | M | € |
| B11 | H-style harness CLASSIC RPRK010 | (XS) S | M | M | M | € |
| B12 | Immo. Waistcoat TBX RPRK039 | (XS) S | M | M | M | € |
| B13 | H-style harness TBX RPRK040 | (XS) S | M | M | M | € |





Pelvis / Abduction



| No. | Desc. / Codes | T4C MINI | T4C STD | T4C MAXI | T4C SM42 | |
|-----|--|----------|---------|----------|----------|---|
| C1 | FLEXI-pads set RPRK033UN | UNI | UNI | UNI | UNI | € |
| C2 | FLEXI-pads set with extended thigh guides RPRK04801 | | | UNI | UNI | € |
| C3 | 2-point immobilisation belt CLASSIC RPRK005 | S | M | M | M | € |
| C4 | Pelvis fixation RPRK076 | S | S | M | M | € |
| C5 | Abduction belt CLASSIC RPRK011 | S | M | M | M | € |
| C6 | Abduction block CLASSIC RPRK007 | S | S | M | M | € |



Legs / Arms





| | | | | | | |
|----|---|-----|-----|-----|-----|---|
| D1 | Footplate angle adjustable RPRK08402 | UNI | UNI | | | € |
| D2 | Footplate angle adjustable RPRK08403 | | | UNI | UNI | € |
| D3 | Footrest-SHORT version RPRK04501 | UNI | | | | € |
| D4 | Feet fixation - pair RPRK008UN | UNI | UNI | UNI | UNI | € |

| No. | Desc. / Codes | T4C MINI | T4C STD | T4C MAXI | T4C SM42 | |
|-----|--|----------|---------|----------|----------|---|
| D5 | Extended side protection CLASSIC RPRK04201 | UNI | UNI | UNI | UNI | € |
| D6 | Padded frame cover with side protection RPRK051UN | UNI | UNI | UNI | UNI | € |
| D7 | Work board CLASSIC RPRK01710 | UNI | UNI | UNI | UNI | € |
| D8 | Work board BASIC RPRK08001 | UNI | UNI | UNI | UNI | € |
| D9 | Footplate soft padding RPRK099 | UNI | UNI | UNI | UNI | € |





Protection / transport



| | | | | | | |
|-----|--|---------|---------|---------|---------|---|
| E1 | Grab bar | UNI(SE) | UNI(SE) | UNI(SE) | UNI(SE) | € |
| E2 | Side cushions | UNI(SE) | UNI(SE) | UNI(SE) | UNI(SE) | € |
| E3 | Padded frame cover RPRK01901 | UNI | UNI | UNI | UNI | € |
| E4 | Belt protector W25 RPRK034UN | UNI | UNI | UNI | UNI | € |
| E5 | Belt protector W40 RPRK035UN | UNI | UNI | UNI | UNI | € |
| E6 | Canopy CLASSIC RPRK00107 | UNI | UNI | UNI | UNI | € |
| E7 | Rain cover RPRK01301 | M | M | | | € |
| E8 | Rain cover RPRK01302 | | | L | L | € |
| E9 | Basket (textile) RPRK02102 | UNI | UNI | UNI | | € |
| E10 | Basket (textile) RPRK02106 | | | | UNI | € |
| E11 | Basket-HD (textile) - (FSCW only) RPRK03701 | UNI | UNI | UNI | | € |
| E12 | Bag STYLE Plus RPRK02203 | UNI | UNI | UNI | UNI | € |
| E13 | Footmuff - Summer RPRK014 | S | M | M | L | € |
| E14 | Footmuff - Winter RPRK015 | S | M | M | L | € |



Specialties / interfacing

| No. | Desc. / Codes | T4C MINI | T4C STD | T4C MAXI | T4C SM42 | |
|-----|--|----------|---------|----------|----------|---|
| F1 | FCSW-with swivel lock (set) RPRK41701 | UNI | UNI | UNI | UNI | € |
| F2 | Additionally with rear AIR wheels RPRK40208 | UNI | UNI | UNI | UNI | € |
| F3 | Front big fixed wheels | UNI | UNI | UNI | UNI | € |
| F4 | Additionally with all AIR wheels RPRK40209 | UNI | UNI | UNI | UNI | € |
| F5 | High-low indor frame CLASSIC RKC101YYY | UNI | UNI | UNI | UNI | € |
| F6 | Oxygen tank holder - vertical 6,5 RPRK09711 | UNI | UNI | UNI | UNI | € |
| F7 | IV pole holder RPRK09811 | UNI | UNI | UNI | UNI | € |
| F8 | Drum-brake adapt.(w.hand control) RPRK41501 | | | UNI | | € |
| F9 | Drum-brake adapt.(w.hand control) RPRK41502 | | | | UNI | € |





Sizes:

- XS - extra small
- S - small
- M - medium
- L - large
- XL - extra large
- UNI - universal size (one size only)

Other notes:

- (SE) Standard Equipment
- (empty position) - accessory is not available