

ORDER FORM

TAMPA




patron[®]
creative in care

Release date 2021/04

Company	Invoice address	Order number / reference	
Street		Customer reference	
E-Mail		Date	
City, Zip code		Offer	
Telephone		Order	
Company	Delivery address		
Street	City, Zip code		

Standard equipment:	Size / order codes		
	STD	MAXI	SM 42
<ul style="list-style-type: none"> • seat depth adjustment • height adjustable backrest • width adjustable seat by FlexiPads set (optimal equipment) • seat height and tilt in space function • angle and height adjustable footrest • height adjustable armrest • rear castor wheels with brake • telescopic wheeled base with length adjustment • grab bar, armrest, push handle. 	TMP101YYY	TMP101XYY	TMP101S2Y
	€	€	€

Backrest height [mm]	520-650	570-680	570-680
Seat depth [mm]	280-310	340-430	340-430
Seat width [mm]	270/340	300/370	350/420
Seat to floor [mm]	530-580	530-580	530-580
Seat to footplate [mm]	260-400	260-400	260-500
Backrest angle	0° - 180°	0° - 160°	0° - 135°
Seat angle	+/-27°	+/-25°	+/-25°
Maximum load	35,0 kg	50,0 kg	75,0 kg

colour frame / pattern colour	
frame	K2 - Anthracite
fabric patterns	R43  Dark lime green
	R44  Levander blue
	R45  Ruby red

TAMPA Accessories

Positioning and Fixation

Head

No.	Desc. / Codes	TMP STD	TMP MAXI	TMP SM42	
A1	Headrest CLASSIC RPRK002	M	M	M	€
A2	Headwings CLASSIC RPRK049	M	M		€
A3	Neckrest CLASSIC RPRK028	M	M	M	€
A4	Headwings RPRK091			UNI	€
A5	Backrest ext.for Headwigs RPRK09204			UNI	€
A6	Headrest contoured, height+angle adj. RPRK093			UNI	€





Trunk

B1	Side trunkrest CLASSIC RPRK003	M	M	L	XL	€
B2	Chest immo. belt CLASSIC RPRK012	M	M	L		€
B3	Hip rest CLASSIC RPRK006	M	L	L		€
B4	Lumbar support W25 RPRK04302	UNI				€
B5	Lumbar support W35 RPRK04303		UNI	UNI		€
B6	5-point immobilisation belt CLASSIC RPRK047	M	M	L		€
B7	4-point immobilisation belt CLASSIC RPRK004	M	M	L		€
B8	Chest-belt addition to 4P-belt RPRK041	UNI	UNI	UNI		€
B9	Immo. Waistcoat CLASSIC RPRK009	M	M	M		€
B10	H-style harness CLASSIC RPRK010	M	M	M		€
B11	Immo. Waistcoat TBX RPRK039	M	M	M		€
B12	H-style harness TBX RPRK040	M	M	M		€








Pelvis / Abduction

No.	Desc. / Codes	TMP STD	TMP MAXI	TMP SM42	
C1	FLEXI-pads set RPRK033UN	UNI	UNI	UNI	€
C2	FLEXI-pads set with extended thigh guides RPRK04801		UNI	UNI	€
C3	2-point immobilisation belt CLASSIC RPRK005	M	M	M	€
C4	Pelvis fixation RPRK076	S	M	M	€
C5	Abduction belt CLASSIC RPRK011	M	M	M	€
C6	Abduction block CLASSIC RPRK007	S	M	M	€




Legs / Arms

No.	Desc. / Codes	TMP STD	TMP MAXI	TMP SM42	
D1	Footplate angle adjustable RPRK08402	UNI			€
D2	Footplate angle adjustable RPRK08403		UNI	UNI	€
D3	Feet fixation - pair RPRK008UN	UNI	UNI	UNI	€
D4	Footrest SHORT version RPRK04501	UNI			€
D5	Footplate soft padding RPRK099	M	L	L	€
D6	Calf rest (textile) RPRK103	UNI	UNI	UNI	€
D7	Work board CLASSIC RPRK01710	UNI	UNI	UNI	€
D8	Work board BASIC RPRK08001	UNI	UNI	UNI	€



Protection / transport

E1	Grab bar	UNI (SE)	UNI (SE)	UNI (SE)	€
E2	Side cushions	UNI (SE)	UNI (SE)	UNI (SE)	€
E3	Belt protector W25 RPRK034	UNI	UNI	UNI	€
E4	Belt protector W40 RPRK035	UNI	UNI	UNI	€



Sizes:

XS - extra small
S - small
M - medium
L - large
XL - extra large
UNI - universal size (one size only)

Other notes:

(SE) Standard Equipment
(empty position) - accessory is not available