

Please contact Invacare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@invacare.com
 Or visit our website : www.invacare.eu.com



Küschall K-Series



Active, lightweight aluminum wheelchair, rigid frame, folding and angle-adjustable backrest



- Description:**
- Aluminum frame with hydroformed tubes
 - Customization kit with premium decals
 - Backrest adjustable in angle and folding on the seat with locking as standard
 - Max user weight: 130 kg
 - Complete weight: from 7.8 kg.
 - Transport weight: from 5.1 kg.
 - Warranty: 2 years excluding wearing parts

Discover our new tool to help calculate overall dimensions and find the visual configurator to make your configuration!

ALL FRAME COLORS & DECALS/STYLINGS FOR FREE UNTIL END OF JUNE!

Date :

Purchase Order :

Customer Code :

Distributor Stamp :

We reserve the right to deliver a + or - tolerance of +/- 10 mm on all measurements We reserve the right to make technical changes without notice


Please make your selection with a cross in the box

◆ Standard feature


○ Free of charge option





CODE	DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros	
DKA0067	<input type="checkbox"/> küschall® K-Series®			
SPECIAL CONSTRUCTION				
DKA0000	<input type="checkbox"/> Myküschall special construction		On Request	
PACKS POSITIONNEMENT ET PRÉVENTION D'ESCARRES				
DKA0220	<input type="checkbox"/> Matrx MX2 & Libra positioning and pressure ulcer prevention pack including: - Matrx MX2 (size of your choice DKA2300 to DKA2315) - Matrx Libra cushion (DKA2720) - Bodypoint® Evoflex® 2-point hip belt (size of your choice DKA7810 to DKA7830)	NC: Seat width 360 mm possible only with depth 400/425 mm, 1310, 1315, 1350, 1355, 1410, 1420	1222,23	Save 15% by selecting this PACK, i.e. €277,53
DKA0230	<input type="checkbox"/> Matrx Elite E2 & Libra positioning and pressure ulcer prevention pack including: - Matrx Elite E2 / Elite Deep E2 backrest (size of your choice DKA2320 to DKA2350) - Matrx Libra cushion (DKA2720) - Bodypoint® Evoflex® 2-point hip belt (size of your choice DKA7810 to DKA7830)	NC: Seat width 360 mm possible only with depth 400/425 mm, 1310, 1315, 1410, 1420	910,24	Save 10% by selecting this PACK, i.e. €106,70
SEATING CHOICE				
FRAME				
DKA1040	<input type="checkbox"/> Frame Angle 75° To The Ground Aluminium	NC : 1475, 1490, 1495, 7240		◆
DKA1050	<input type="checkbox"/> Frame Angle 75° To The Ground Aluminium - Tapered Front Frame	NC : 1475, 1490, 1495, 1740, 1745, 1750, 3230, 7240		○
DKA1060	<input type="checkbox"/> Frame Angle 75° To The Ground Aluminium +50Mm	NC : 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760	145,97	<p>Châssis 90° : 1070</p>
DKA1070	<input type="checkbox"/> Frame Angle 90° To The Ground Aluminium	NC : 1475, 1490, 1495, 7240	217,06	
DKA1080	<input type="checkbox"/> Frame Angle 90° To The Ground Aluminium - Tapered Front Frame	NC : 1475, 1490, 1495, 1740, 1745, 1750, 3230, 5300, 5350, 5400, 5450, 5500, 7240	217,06	
DKA1090	<input type="checkbox"/> Frame Angle 90° To The Ground Aluminium +50 Mm	NC : 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760	217,06	

SUSPENSION				
DKA1200	<input type="checkbox"/>	40 - 65 kg	NC : 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430, 6950	336,49
DKA1210	<input type="checkbox"/>	65 - 85 kg	NC : 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430, 6950	336,49
DKA1220	<input type="checkbox"/>	85 - 100 kg	NC : 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430, 6950	336,49
SEAT WIDTH (SW)				
DKA1310	<input type="checkbox"/>	320 mm	NC : 2000, 2050, 2300 to 2357, 2610, 2630, 2710, 2720, 2740, 2742, 2744, 2746, 2748, 3405, 3420	○
DKA1315	<input type="checkbox"/>	340 mm	NC : 2300 to 2357, 2710, 2720, 2740, 2742, 2744, 2746, 2748, 3405, 3420	○
DKA1320	<input type="checkbox"/>	360 mm	NC : 2364, 2366, 2368, 2372, 2710, 2742, 2744, 2746, 2748	○
DKA1325	<input type="checkbox"/>	380 mm	NC : 2366, 2368, 2370, 2372	○
DKA1330	<input type="checkbox"/>	400 mm	NC : 2360, 2366, 2368, 2370, 2372	○
DKA1335	<input type="checkbox"/>	420 mm	NC : 2360, 2362, 2368, 2370, 2372	○
DKA1340	<input type="checkbox"/>	440 mm	NC : 2210, 2360, 2362, 2364, 2370, 2372, 2630, 2740, 2742	○
DKA1345	<input type="checkbox"/>	460 mm	NC : 2210, 2360, 2362, 2364, 2372, 2630, 2740, 2742, 2744	○
DKA1350	<input type="checkbox"/>	480 mm	NC : 2210, 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2630, 2740, 2742, 2744, 3230	○
DKA1355	<input type="checkbox"/>	500 mm	NC : 2210, 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2368, 2630, 2740, 2742, 2744, 2746, 3210, 3230	○
SEAT DEPTH (SD)				
DKA1410	<input type="checkbox"/>	350 mm	NC : 1500, 1510, 2630, 2710, 2720, 2740 to 2746, 2748, 2860, 7810, 7820, 7830	○
DKA1420	<input type="checkbox"/>	375 mm	NC : 2710, 2720, 2742, 2744, 2746, 2748, 7810, 7820, 7830	○
DKA1435	<input type="checkbox"/>	400 mm		○
DKA1450	<input type="checkbox"/>	425 mm	NC : 2740, 2746	○
DKA1465	<input type="checkbox"/>	450 mm	NC : 2630, 2740, 2742	○
DKA1475	<input type="checkbox"/>	475 mm	NC : 2630, 2740, 2742, 2744, 2748	○
DKA1490	<input type="checkbox"/>	500 mm	NC : 2630, 2710, 2740, 2742, 2744, 2746, 2748	○
DKA1495	<input type="checkbox"/>	525 mm	NC : 2610, 2630, 2710, 2740, 2742, 2744, 2746, 2748	○
SEAT POSITION (SP)				
DKA1500	<input type="checkbox"/>	Position 1 = 140 mm (MOST active)		○
DKA1510	<input type="checkbox"/>	Position 2 = 125 mm		○
DKA1520	<input type="checkbox"/>	Position 3 = 110 mm		◆
DKA1530	<input type="checkbox"/>	Position 4 = 95 mm		○
DKA1540	<input type="checkbox"/>	Position 5 = 80 mm (Least stable)		○













BACK HEIGHT (BH)						
DKA1600	<input type="checkbox"/>	Backrest height according to chosen Matrix backrest		NC : 2100, 2210, 2220, 2230	○	
DKA1605	<input type="checkbox"/>	270 mm		NC : 2150, 2240, 2430, 0703, 7600	○	
DKA1610	<input type="checkbox"/>	285 mm		NC : 2150, 2240, 7600	○	
DKA1615	<input type="checkbox"/>	300 mm	DKA1620 <input type="checkbox"/>	315 mm	NC : 2150, 2240, 7600	○
DKA1625	<input type="checkbox"/>	330 mm	DKA1630 <input type="checkbox"/>	345 mm	NC : 2150, 2240, 7600	○
DKA1635	<input type="checkbox"/>	360 mm		NC : 2150, 2210, 2240, 7600	○	
DKA1640	<input type="checkbox"/>	375 mm		NC : 2150, 2210, 2240, 7600	◆	
DKA1645	<input type="checkbox"/>	390 mm		NC : 2150, 2210, 2240, 7600	○	
DKA1650	<input type="checkbox"/>	405 mm	DKA1655 <input type="checkbox"/>	420 mm	NC : 2150, 2210, 2240	○
DKA1660	<input type="checkbox"/>	435 mm	DKA1665 <input type="checkbox"/>	450 mm	NC : 2150, 2210, 2240	○
DKA1670	<input type="checkbox"/>	465 mm	DKA1675 <input type="checkbox"/>	480 mm	NC : 2150, 2210, 2240	○
Front seat height (FSTF) - see table 1 of front wheel compatibility and front seat height						
DKA1740	<input type="checkbox"/>	450 mm		NC : 5300, 5350, 5400, 5450, 5500	○	
DKA1745	<input type="checkbox"/>	460 mm		NC : 5300, 5350, 5400, 5450, 5500	○	
DKA1750	<input type="checkbox"/>	470 mm			○	
DKA1755	<input type="checkbox"/>	480 mm			○	
DKA1760	<input type="checkbox"/>	490 mm			○	
DKA1765	<input type="checkbox"/>	500 mm			◆	
DKA1770	<input type="checkbox"/>	510 mm			○	
DKA1775	<input type="checkbox"/>	520 mm			○	
REAR SEAT-TO-FLOOR (RSTF)						
DKA1815	<input type="checkbox"/>	390 mm		NC : 3490, 3495, 6030, 6040, 6420, 6430	○	
DKA1820	<input type="checkbox"/>	400 mm		NC : 3490, 3495, 6030, 6040, 6430	○	
DKA1825	<input type="checkbox"/>	410 mm		NC : 3490, 3495	○	
DKA1830	<input type="checkbox"/>	420 mm			○	
DKA1835	<input type="checkbox"/>	430 mm			◆	
DKA1840	<input type="checkbox"/>	440 mm	DKA1845 <input type="checkbox"/>	450 mm	○	
DKA1850	<input type="checkbox"/>	460 mm	DKA1855 <input type="checkbox"/>	470 mm	○	
DKA1860	<input type="checkbox"/>	480 mm			○	
DKA1865	<input type="checkbox"/>	490 mm		NC : 6410	○	
LOWER LEG LENGTH (LLL) - FOR LLL < 390MM AND DEPENDING ON CONFIGURATION, FOOTPLATE MAY BE UP MOUNTED, SEE TABLE 2						
DKA1918	<input type="checkbox"/>	290 mm			○	
DKA1920	<input type="checkbox"/>	300 mm	DKA1922 <input type="checkbox"/>	310 mm	○	
DKA1924	<input type="checkbox"/>	320 mm	DKA1926 <input type="checkbox"/>	330 mm	○	
DKA1928	<input type="checkbox"/>	340 mm	DKA1930 <input type="checkbox"/>	350 mm	○	
DKA1932	<input type="checkbox"/>	360 mm	DKA1934 <input type="checkbox"/>	370 mm	○	
DKA1936	<input type="checkbox"/>	380 mm	DKA1938 <input type="checkbox"/>	390 mm	○	
DKA1940	<input type="checkbox"/>	400 mm	DKA1942 <input type="checkbox"/>	410 mm	○	
DKA1944	<input type="checkbox"/>	420 mm			○	
DKA1946	<input type="checkbox"/>	430 mm			◆	
DKA1948	<input type="checkbox"/>	440 mm	DKA1950 <input type="checkbox"/>	450 mm	○	
DKA1952	<input type="checkbox"/>	460 mm	DKA1954 <input type="checkbox"/>	470 mm	○	
DKA1956	<input type="checkbox"/>	480 mm			○	
BACKREST						
BACKREST MATERIAL						
DKA2000	<input type="checkbox"/>	Aluminium		NC : SW 320 mm	◆	
DKA2050	<input type="checkbox"/>	Titanium = - 80 g (with aluminum push handles)		NC : SW 320 mm, 2150, 2430, 7600	345,02	
BACKREST TYPE						
DKA2100	<input type="checkbox"/>	Standard		NC : 2240, 2300 to 2357, 2440	◆	
DKA2150	<input type="checkbox"/>	BACKREST WITH LOWER BACKBRACE, BACKREST TUBES HEIGHT IS 285 MM		NC : 2210, 2220, 2230, 2430	○	



BACK UPHOLSTERY					
DKA2210	<input type="checkbox"/>	LIGHT BACK UPHOLSTERY	NC : 2300 to 2357, 2430, 2440 and SW > 420 mm / BH > 345 mm	178,20	 2230
DKA2220	<input type="checkbox"/>	TENSION ADJUSTABLE PREMIUM BLACK BACK UPHOLSTERY WITH LOOP BANDS	NC : 2300 to 2357, 2440	◆	
DKA2230	<input type="checkbox"/>	AERO BACK, TENSION ADJUSTABLE	NC : 2300 to 2357, 2440	258,77	
DKA2240	<input type="checkbox"/>	WITHOUT BACKREST UPHOLSTERY, SELECT IF YOU WANT TO MOUNT POSITIONNING BACKRESTS SUCH AS MATRX		○	
MATRX BACKRESTS - HEIGHT REFERS TO ACTUAL HEIGHT OF THE MATRX BACKREST, NOT POSSIBLE WITH SW ≤ 340 MM					
DKA2300	<input type="checkbox"/>	Matrx MX2, Height 150 Mm With Quick Release Hardware	NC : 2370, 2372	927,01	 2300, 2305, 2310, 2315
DKA2305	<input type="checkbox"/>	Matrx MX2, Height 230 Mm With Quick Release Hardware	NC : 2370, 2372	927,01	
DKA2310	<input type="checkbox"/>	Matrx MX2, Height 305 Mm With Quick Release Hardware	NC : 2370, 2372	927,01	
DKA2315	<input type="checkbox"/>	Matrx MX2, Height 405 Mm With Quick Release Hardware	NC : 2370, 2372	927,01	
DKA2316	<input type="checkbox"/>	Elite Deep E2 , Height 160 Mm With Mini-Set Mounting Hardware		500,47	
DKA2317	<input type="checkbox"/>	Elite Deep E2 , Height 210 Mm With Mini-Set Mounting Hardware		500,47	<div style="border: 1px solid blue; padding: 5px; color: blue; font-size: small;"> MX2: only 140g more than a tension-adjustable backrest! </div>
DKA2318	<input type="checkbox"/>	Elite Deep E2 , Height 260 Mm With Mini-Set Mounting Hardware		500,47	
DKA2320	<input type="checkbox"/>	Elite Deep E2 , Height 310 Mm With Mini-Set Mounting Hardware		500,47	
DKA2322	<input type="checkbox"/>	Elite Deep E2, Height 360 Mm With Mini-Set Mounting Hardware		500,47	
DKA2325	<input type="checkbox"/>	Elite Deep E2, Height 410 Mm With Mini-Set Mounting Hardware		500,47	
DKA2330	<input type="checkbox"/>	Elite Deep E2, Height 460 Mm With Easy-Set Mounting Hardware		500,47	
DKA2331	<input type="checkbox"/>	Elite Deep E2, Height 510 Mm With Easy-Set Mounting Hardware		500,47	
DKA2332	<input type="checkbox"/>	Elite E2, Height 160 Mm With Mini-Set Mounting Hardware		500,47	
DKA2334	<input type="checkbox"/>	Elite E2, Height 210 Mm With Mini-Set Mounting Hardware		500,47	
DKA2335	<input type="checkbox"/>	Elite E2, Height 260 Mm With Mini-Set Mounting Hardware		500,47	
DKA2340	<input type="checkbox"/>	Elite E2, Height 310 Mm With Mini-Set Mounting Hardware		500,47	 2316 to 2331
DKA2345	<input type="checkbox"/>	Elite E2, Height 360 Mm With Mini-Set Mounting Hardware		500,47	
DKA2350	<input type="checkbox"/>	Elite E2, Height 410 Mm With Mini-Set Mounting Hardware		500,47	
DKA2355	<input type="checkbox"/>	Elite E2, Height 460 Mm With Easy-Set Mounting Hardware		500,47	
DKA2357	<input type="checkbox"/>	Elite E2, Height 510 Mm With Easy-Set Mounting Hardware		500,47	
MATRX BACKRESTS WIDTH CHOICE - See Table 3					
DKA2360	<input type="checkbox"/>	360 mm - indicative bust 265 mm		○	 2332 to 2357
DKA2362	<input type="checkbox"/>	380 mm - indicative bust 290 mm		○	
DKA2364	<input type="checkbox"/>	410 mm - indicative bust 315 mm		○	
DKA2366	<input type="checkbox"/>	430 mm - indicative bust 340 mm		○	
DKA2368	<input type="checkbox"/>	460 mm - indicative bust 370 mm		○	
DKA2370	<input type="checkbox"/>	480 mm - indicative bust 390 mm		○	
DKA2372	<input type="checkbox"/>	510 mm - indicative bust 420 mm		○	

PUSH HANDLES				
DKA2400	<input type="checkbox"/>	WITHOUT PUSH HANDLES		◆
DKA2410	<input type="checkbox"/>	PUSH HANDLES FIXED, SHORT		○
DKA2420	<input type="checkbox"/>	PUSH HANDLES FIXED, LONG		○
DKA2430	<input type="checkbox"/>	HEIGHT ADJUSTABLE, REAR SET	NC : BH 285 mm to 330 mm + 3480, 3480	212,32
DKA2440	<input type="checkbox"/>	HEIGHT ADJUSTABLE, INTEGRATED		254,03
DKA2450	<input type="checkbox"/>	FOLDABLE PUSH HANDLES	NC : BH 315 mm to 330 mm + 3480	213,27
BACKREST ANGLE (BA)				
DKA2520	<input type="checkbox"/>	78°		○
DKA2540	<input type="checkbox"/>	82°		◆
DKA2560	<input type="checkbox"/>	86°		○
DKA2580	<input type="checkbox"/>	90°	DKA2590 <input type="checkbox"/>	94°
				○
SEAT UPHOLSTERY				
DKA2600	<input type="checkbox"/>	STANDARD SEAT UPHOLSTERY - CENTRAL LOOP BAND		◆
DKA2610	<input type="checkbox"/>	TENSION ADJUSTABLE WITH ONE CENTRAL LOOP BAND		92,89
DKA2630	<input type="checkbox"/>	LIGHT SEAT UPHOLSTERY		66,35
SEATING				
SEAT CUSHION				
DKA2700	<input type="checkbox"/>	STANDARD SEAT CUSHION,60 MM, DOUBLE FOAM WITH VISCO TOP	NC : 2810, 2820, 2830, 2840, 2850	101,42
DKA2710	<input type="checkbox"/>	FLO-TECH LITE	NC : 2810, 2820, 2830, 2840, 2850	118,48
DKA2720	<input type="checkbox"/>	MATRIX LIBRA CUSHION, BLACK SMOOTH COVER	NC : largeur 360 mm possible uniquement avec profondeur 400 and 425 mm	293,84
DKA2740	<input type="checkbox"/>	CLASSIC VISCO CUSHION 375 X 400 MM	NC : 2810, 2820, 2830, 2840, 2850	76,78
DKA2742	<input type="checkbox"/>	CLASSIC VISCO CUSHION 405 X 400 MM	NC : 2810, 2820, 2830, 2840, 2850	65,40
DKA2744	<input type="checkbox"/>	CLASSIC VISCO CUSHION 405 X 400 MM	NC : 2810, 2820, 2830, 2840, 2850	65,40
DKA2746	<input type="checkbox"/>	CLASSIC VISCO CUSHION 455 X 400 MM	NC : 2810, 2820, 2830, 2840, 2850	76,78
DKA2748	<input type="checkbox"/>	CONTOUR VISCO NG CUSHION	NC : 2810, 2820, 2830, 2840, 2850	76,78
SEAT BAG				
DKA2800	<input type="checkbox"/>	SEAT BAG, CLICK-IN AND REMOVABLE		81,52
MATRIX LIBRA CUSHION ACCESSORIES				
DKA2810	<input type="checkbox"/>	FLUID SUPPLEMENT SAC (PAIR)		91,00
DKA2820	<input type="checkbox"/>	TAPERED WEDGE PAIR		41,71
DKA2830	<input type="checkbox"/>	ABDUCTOR		47,39
DKA2840	<input type="checkbox"/>	LATERAL SIDE WEDGE		41,71
DKA2850	<input type="checkbox"/>	ANTERIOR FRONT WEDGE		54,98
Base rigide				
DKA2860	<input type="checkbox"/>	SEAT RIGIDIZER		56,87
SEAT UPHOLSTERY BAGS				
DKA2910	<input type="checkbox"/>	BIG BAG - TO FIX UNDER THE SEAT		61,61



FOOTREST					
FOOTREST TYPE (TUBE / PLATE)					
DKA3200	<input type="checkbox"/>	FOOTREST TUBE, RIGID - ALUMINIUM NEW PROFILE		◆	 3200
DKA3210	<input type="checkbox"/>	FOOTREST TUBE, RIGID - TITANIUM		246,45	
DKA3220	<input type="checkbox"/>	CARBON FOOTPLATE, ANGLE ADJUSTABLE	NC : 3290	264,45	
DKA3230	<input type="checkbox"/>	FOOTPLATE FLIP-UP TO SIDE, ANGLE ADJUSTABLE NOT POSSIBLE WITH TAPERED FRAME	NC : 3290	458,77	
FOOTREST COVER					
DKA3290	<input type="checkbox"/>	FOOTREST COVER			 3305
DKA3290	<input type="checkbox"/>	FOOTREST COVER		113,74	
CALF STRAP MOUNTED ON THE FRAME					
DKA3300	<input type="checkbox"/>	CALF STRAP MOUNTED ON THE FRAME		◆	3305
DKA3305	<input type="checkbox"/>	CALF STRAP WITH BRACKET SUPPORT		76,78	
CLOTHES GUARDS / MUDGUARDS					
CLOTHES GUARDS					
DKA3405	<input type="checkbox"/>	ECONOMY CLOTHES GUARDS, FLIP-BACK - NEW PART		◆	 3440
DKA3420	<input type="checkbox"/>	CARBON CLOTHES GUARDS, FLIP BACK - NEW DESIGN		167,77	
MUDGUARDS					
DKA3440	<input type="checkbox"/>	CARBON MUDGUARDS, REMOVABLE- NEW DESIGN	NC : 3490, 3495		 3480
DKA3450	<input type="checkbox"/>	ECONOMY MUDGUARDS, REMOVABLE - NEW DESIGN	NC : 3490, 3495	○	
ARMREST / SIDEREST					
DKA3480	<input type="checkbox"/>	TUBULAR ARMREST, NON-LOCKING, HEIGHT ADJUSTABLE, REMOVABLE, SWING-AWAY, COMFORT PADDED	NC : BH < 315 mm + 2450, BH <315 + 2430, 6950		 3490
DKA3490	<input type="checkbox"/>	SIDEREST, HEIGHT ADJUSTABLE, SHORT, NON-LOCKING	NOT POSSIBLE WITH RSTF < 420 MM	253,08	
DKA3495	<input type="checkbox"/>	SIDEREST, HEIGHT ADJUSTABLE, LONG, NON-LOCKING	NOT POSSIBLE WITH RSTF < 420 MM	277,73	
FRAME					
FRAME COLOR					
DKA4000	<input type="checkbox"/>	SPACE BLACK (RAL 9005)		○	 EXTRA LEAD TIME +10 DAYS
DKA4010	<input type="checkbox"/>	ANTHRACITE GREY (RAL 7016)		○	
DKA4020	<input type="checkbox"/>	PURE WHITE (RAL 9016)	NC : 4220	○	
DKA4030	<input type="checkbox"/>	CARMINE RED (RAL 3002)		○	
DKA4040	<input type="checkbox"/>	AZURE BLUE (RAL 5017)		○	
DKA4060	<input type="checkbox"/>	FLASHY YELLOW (RAL 1016)		FOR FREE!	
DKA4070	<input type="checkbox"/>	ATOMIC VIOLET - specific color mix, no RAL definition			
DKA4100	<input type="checkbox"/>	BLACK GLOSSY (RAL 9005)			
DKA4130	<input type="checkbox"/>	SILVER GREY (RAL 9006)			
DKA4140	<input type="checkbox"/>	ROYAL BLUE - specific color mix, no RAL definition			
DKA4150	<input type="checkbox"/>	RUBY RED (RAL 3003)			
DKA4170	<input type="checkbox"/>	DYNAMITE ORANGE - specific color mix, no RAL definition			
DKA4180	<input type="checkbox"/>	FOREST GREEN (RAL 6029)			
DKA4050	<input type="checkbox"/>	Special colors on demand only : - Manufacturing time + 10 days: for a color already offered on the other ranges of Invacare manual wheelchairs. WHERE - Manufacturing time + 15 days: for any other request (subject to feasibility) - specify RAL number		473,93	
STYLING					
DKA4190	<input type="checkbox"/>	BOLD	NC : 4300, 4310	FOR FREE!	 BOLD  TAG
DKA4195	<input type="checkbox"/>	TAG	NC : 4200 to 4250	◆	
BOLD COLOR					
DKA4200	<input type="checkbox"/>	BLACK		○	
DKA4210	<input type="checkbox"/>	DARK GREY		○	
DKA4220	<input type="checkbox"/>	LIGHT GREY		○	
DKA4230	<input type="checkbox"/>	RED		○	
DKA4240	<input type="checkbox"/>	BLUE		○	
DKA4250	<input type="checkbox"/>	GOLD		○	

TAG COLOR			
DKA4300	<input type="checkbox"/>	BLACK	○
DKA4310	<input type="checkbox"/>	RED	○
FRONT WHEELS			
DKA5000	<input type="checkbox"/>	3" SUPERSPORT	○
DKA5050	<input type="checkbox"/>	4" EVERYDAY WHEEL WITH BLACK TYRE	◆
DKA5100	<input type="checkbox"/>	4" COMFORT WHEEL WITH BLACK TYRE	○
DKA5150	<input type="checkbox"/>	4" STAREC BLACK HUB, BLACK TYRE	75,83
DKA5200	<input type="checkbox"/>	4" FROGLEGS FRONT WHEEL WITH BLACK HUB, BLACK TYRE	121,33
DKA5250	<input type="checkbox"/>	4" FROGLEGS FRONT WHEEL WITH SILVER HUB, BLACK TYRE	121,33
DKA5300	<input type="checkbox"/>	5" EVERYDAY WITH BLACK RUBBER HUB, BLACK TYRE	○
DKA5350	<input type="checkbox"/>	5" COMFORT WHEEL BLACK RUBBER HUB, BLACK TYRE	○
DKA5400	<input type="checkbox"/>	5" STAREC BLACK HUB, BLACK TYRE	103,32
DKA5450	<input type="checkbox"/>	5" FROGLEGS FRONT WHEEL WITH BLACK HUB, BLACK TYRE	147,87
DKA5500	<input type="checkbox"/>	5" FROGLEGS FRONT WHEEL WITH SILVER HUB, BLACK TYRE	147,87
REAR WHEEL FIXATION			
DKA6020	<input type="checkbox"/>	STANDARD REAR WHEEL POSITION	NC : 6100 to 6140, 6230, 6240
DKA6030	<input type="checkbox"/>	LIGHT REAR WHEEL EXTENSION (POSITION 1 TO 3)	NC : 6070, 6130, 6140, 6230, 6240
DKA6040	<input type="checkbox"/>	REAR WHEEL EXTENSION (5 POSITIONS)	NC : 6070
CAMBER AXLE TYPE - 1st AXLE			
DKA6050	<input type="checkbox"/>	VARIO-AX, ADJUSTABLE CAMBER 1 / 3 / 7 / 10°	NC : 6200, 6210, 6220, 6270, 6950, 7230
DKA6060	<input type="checkbox"/>	STANDARD AXLE	◆
DKA6070	<input type="checkbox"/>	CARBON AXLE	NC : 6950
POSITION OF 1st AXLE			
DKA6100	<input type="checkbox"/>	1 (compulsory for the use of a 2nd axle)	○
DKA6110	<input type="checkbox"/>	2	NC : 6230, 6240
DKA6120	<input type="checkbox"/>	3	NC : 6230, 6240
DKA6130	<input type="checkbox"/>	4	NC : 6230, 6240
DKA6140	<input type="checkbox"/>	5	NC : 6230, 6240
CAMBER ANGLE 1st AXLE			
DKA6200	<input type="checkbox"/>	0°	NC : 6270
DKA6210	<input type="checkbox"/>	3°	◆
DKA6220	<input type="checkbox"/>	6°	NC : 6240, 6250, 6950
SECOND CAMBER AXLE TYPE - SECOND AXLE ONLY IN THE 5TH POSITION FROM THE FRONT			
DKA6230	<input type="checkbox"/>	STANDARD	NC : 6950
DKA6240	<input type="checkbox"/>	VARIO-AX, ADJUSTABLE CAMBER 1 / 3 / 7 / 10°	NC : 6250, 6260, 6270, 6950, 7230
CAMBER ANGLE 2nd AXLE - MAX. 3° CAMBER DIFFERENCE WITH 1ST AXLE ; IF ONLY ONE OF THE AXLES IS A VARIO-AX, 2ND AXLE ONLY IN 0 OR 3°			
DKA6250	<input type="checkbox"/>	0°	○
DKA6260	<input type="checkbox"/>	3°	○
DKA6270	<input type="checkbox"/>	6°	○
QUAD RELEASE AXLE			
DKA6300	<input type="checkbox"/>	TETRA-CLIP	NC : 6550, 6950
REAR WHEELS SIZE (see the possible diameters per wheel in Type of rear wheels)			
DKA6410	<input type="checkbox"/>	24"	◆
DKA6420	<input type="checkbox"/>	25"	NC : 6500, 6510, 6600, 6610, 6630, 6660, 6680, 6700, 6710, 6720, 6760, 6770, 6900, 6910, 6950
DKA6430	<input type="checkbox"/>	26"	NC : 6500, 6520, 6550, 6600, 6610, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6950, 7240



5000



5200 / 5250 / 5450 / 5500

5350



5050 / 5300

Pos.1-3



Pos.1-5



6300

REAR WHEELS TYPE					
DKA6500	<input type="checkbox"/>	24" ACTIVE WHEELS (BLACK RIM, HUB AND RADIAL SPOKING INCL. ALU HANDRIMS & LIGHT TYRE GREY)	24"	NC : All handrims, 6680, 6900, 6910	◆
DKA6510	<input type="checkbox"/>	KUSCHALL STAREC HUB (BLACK RIM AND HUB, BLACK CROSS SPOKING INCL. ALU HANDRIMS & LIGHT TYRE GREY)	24", 26"		○
DKA6520	<input type="checkbox"/>	HIGH PERFORMANCE (BLACK RIM, BLACK STAREC HUB, BLACK CROSS SPOKING INCL. ALU HANDRIMS & LIGHT TYRE GREY)	24", 25"	NC : 6660	533,65
DKA6550	<input type="checkbox"/>	SPINERGY LXX (BLACK RIM, HUB AND RADIAL SPOKING INCL. ALU HANDRIMS & LIGHT TYRE GREY)	24", 25"	NC : 6660, 6910	771,56
REAR WHEEL TYRE					
DKA6600	<input type="checkbox"/>	MAXXIS LITE GUARD BLACK TYRE (puncture protection)	24"		◆
DKA6625	<input type="checkbox"/>	MAXXIS ALPHA PROTECT BLACK TYRE (puncture protection)	24"		72,04
DKA6630	<input type="checkbox"/>	MBL SPEEDLITE (puncture protection)	24"		72,04
DKA6660	<input type="checkbox"/>	SOLID TYRE, BLACK	24"	NC : 6680, 7230	72,04
HANDRIMS					
DKA6700	<input type="checkbox"/>	ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM	24", 26"		142,18
DKA6710	<input type="checkbox"/>	STEEL, DIST. 3 CM ADJUSTABLE 4CM	24", 26"		113,74
DKA6720	<input type="checkbox"/>	TITANIUM, DIST. 3 CM ADJUSTABLE 4CM	24"		277,73
DKA6740	<input type="checkbox"/>	CARBOLIFE GEKKO , WITH ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM	24", 25"		373,46
DKA6750	<input type="checkbox"/>	CARBOLIFE GEKKO S , WITH ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM	24", 25"		373,46
DKA6760	<input type="checkbox"/>	CARBOLIFE CURVE L WITH ERGONOMIC PROFILE DIST. 3 CM ADJUSTABLE 4CM	24"		282,46
DKA6770	<input type="checkbox"/>	CARBOLIFE QUADRO , WITH ERGONOMIC PROFILE DIST. 3 CM ADJUSTABLE 4CM	24"		373,46
DKA6780	<input type="checkbox"/>	LIGHT ALU, GREY, DIST. 3 CM ADJUSTABLE 4CM	24", 25", 26"		◆
DKA6790	<input type="checkbox"/>	LIGHT ALU, BLACK, DIST. 3 CM ADJUSTABLE 4CM	24", 25", 26"		○



6790 6760



6700 6710



6720



6740


SPOKE GUARDS				
DKA6900	<input type="checkbox"/>	Transparent	24", 26"	200,95
DKA6910	<input type="checkbox"/>	Black	24", 26"	200,95
ALBER ADD-ON KIT - 2 ALBER ANTI-TIPPERS INCLUDED / ONLY FOR 24" REAR WHEELS ALBER ADDS ON KIT + 3° CAMBER AXLE : MAX USER WEIGHT 100 Kg				
DKA6950	<input type="checkbox"/>	E-MOTION & TWION M25 ADAPTER PLATE	NC : 7500, 7510, 7520, 8200, 8210, 8300	579,15
<i>Tip : for perfect compatibility of the brakes with the manual wheels, we advise you to equip them with the same coating as the motorized wheels.</i>				
BRAKES				
DKA7210	<input type="checkbox"/>	PUSH BRAKE, LEVER BENT		◆
DKA7230	<input type="checkbox"/>	PERFORMANCE BRAKE	NC : 7300	76,78
DKA7240	<input type="checkbox"/>	ACTIVE BRAKE	NC : 7300	76,78
DKA7250	<input type="checkbox"/>	LIGHT ACTIVE BRAKE	NC : 7300	145,02
BRAKE EXTENSION LEVER				
DKA7300	<input type="checkbox"/>	EXTENSION LEVER FOR STANDARD BRAKE		76,78
SAFETY OPTIONS				
MAX. 2 ADAPTATIONS PER SIDE (ADAPTATIONS=ANTITIPPER, TIPPER AID, CANE HOLDER, TRANSIT WHEELS)				
ANTI TIPPER				
DKA7500	<input type="checkbox"/>	ANTITIPPER TUBE, RIGHT		130,81
DKA7510	<input type="checkbox"/>	ANTITIPPER TUBE, LEFT		130,81
DKA7520	<input type="checkbox"/>	ANTITIPPER TUBE, BOTH SIDES	NC : 8300	261,61
TRANSPORTATION KIT				
DKA7600	<input type="checkbox"/>	TRANSPORTATION KIT INCL. POSTURE BELT AND LABELS	NC : BH ≤ 400 mm, 7800	57,82
BELT				
DKA7800	<input type="checkbox"/>	STANDARD BELT - PLASTIC BUCKLE		42,65
DKA7810	<input type="checkbox"/>	BODYPOINT POSTURE BELT EVOFLEX, 2 POINTS, SMALL, FLAT-MOUNT WITH PUSH BUTTON		217,06
DKA7820	<input type="checkbox"/>	BODYPOINT POSTURE BELT EVOFLEX, 2 POINTS, MEDIUM, FLAT-MOUNT WITH REHAB LATCH		217,06
DKA7830	<input type="checkbox"/>	BODYPOINT POSTURE BELT EVOFLEX, 2 POINTS, LARGE, FLAT-MOUNT WITH REHAB LATCH		217,06
PASSIVE LIGHTS				
DKA7900	<input type="checkbox"/>	PASSIVE LIGHTS, ASSEMBLED		22,75
ATTENDANT OPTIONS				
FRAME PROTECTION				
DKA8010	<input type="checkbox"/>	FRONT FRAME PROTECTION		44,55
CANE HOLDER				
DKA8200	<input type="checkbox"/>	CANE HOLDER, RIGHT		66,35
DKA8210	<input type="checkbox"/>	CANE HOLDER, LEFT		66,35
TRANSIT WHEELS				
DKA8300	<input type="checkbox"/>	TRANSIT WHEELS		301,42
OTHERS OPTIONS				
DKA8500	<input type="checkbox"/>	AIR PUMP - NEW	NOT POSSIBLE WITH SW 340 MM	33,18
DKA8550	<input type="checkbox"/>	TOOL KIT		22,75
DKA8650	<input type="checkbox"/>	WHEELS BAG		65,40
DKA8680	<input type="checkbox"/>	24" MOUNTAIN BIKE WHEELS WITH BAG -MOUNTAIN BIKE WHEELS AND TIRES WITH ALU GREY HANDRIMS ; QUICK RELEASE AXLES INCLUDED		646,45
DKA8800	<input type="checkbox"/>	BODYPOINT CALF STRAP		104,27
DKA8850	<input type="checkbox"/>	BODYPOINT MOBILITY BAG		104,27



7210 7230 7240



7250



7500 / 7520



Evoflex belt adjustable from 35 mm to 125 mm rearward from the front of the seat sling. 6 positions in 15 mm increments



7810 / 7820 / 7830



8300



8550



8680

RAMPS

1440530	<input type="checkbox"/>	Fiberglass Folding Mono Ramp Folded dimensions: 2000 x 230 mm 16.3 kg - maximum capacity: 300 kg	1 787,68
1654128	<input type="checkbox"/>	Basic ramp 2 telescopic parts in aluminum and non-slip coating Folded dimensions: 1210 x 194 mm 8.8 kg per pair - maximum capacity: 270 kg	402,84
1654129	<input type="checkbox"/>	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - max. : 260 kg	805,69
1654130	<input type="checkbox"/>	Aluminum Folding Mono Ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg	426,54




1440530



1654128



1654129



1654130

Table 1
Frame type/angle/length, rear wheel size, front wheel size

Frame	24" rear wheel			25" rear wheel			26" rear wheel		
	front wheel			front wheel			front wheel		
	3"	4"	5"	3"	4"	5"	3"	4"	5"
75° Short									
75° Standard									
75° +50mm									
90° Short									
90° Standard									
90° +50mm									
75° tapered									
90° tapered									

Table 3
MATRIX BACKREST WIDTH AVAILABLE RANGE RANGE

SEAT WIDTH (SW)	ELITE E2 / ELITE DEEP E2 MATRIX BACKREST WIDTH CHOICE						
	360 MM	380 MM	410 MM	430 MM	460 MM	480 MM	510 MM
320 MM							
340 MM							
360 MM							
380 MM							
400 MM							
420 MM							
440 MM							
460 MM							
480 MM							
500 MM							

SEAT WIDTH (SW)	MX2 MATRIX BACKREST WIDTH CHOICE						
	360 MM	380 MM	410 MM	430 MM	460 MM	480 MM	510 MM
320 MM							
340 MM							
360 MM							
380 MM							
400 MM							
420 MM							
440 MM							
460 MM							
480 MM							
500 MM							

Table 2 :
Frame type/angle/length, rear wheel size, lower leg length, FSTF, footrest type

Frame	24" rear wheel				25" rear wheel				26" rear wheel			
	FSTF	in or bracket	LLL		FSTF	in or bracket	LLL		FSTF	in or bracket	LLL	
			Normal	Flip-up			Normal	Flip-up			Normal	Flip-up
75° Short	430	S	330-420	290-420	460	S	330-430	290-430	470	S	330-440	290-440
	460	S	340-430	300-430	470	S	340-440	300-440	480	S	340-450	300-450
	470	L	350-440	310-440	480	L	350-450	310-450	490	L	350-460	310-460
75° Standard	460	S	370-460	330-460	490	S	370-460	330-460	500	S	370-470	330-470
	490	S	380-460	340-460	500	S	380-470	340-460	510	S	380-480	340-480
	500	L	390-470	350-470	510	L	390-480	350-480	520	L	390-490	350-490
75° +50mm	500	S	390-470	350-470	510	S	390-480	350-480	520	S	390-490	350-490
	510	S	400-480	360-480	520	S	400-490	360-490	530	S	400-500	360-500
	520	L	410-490	370-490	530	L	410-500	370-500	540	L	410-510	370-510
90° Short	430	S	330-410	290-410	460	S	330-420	290-420	470	S	330-430	290-430
	460	S	340-420	300-420	470	S	340-430	300-430	480	S	340-440	300-440
	470	L	350-430	310-430	480	L	350-440	310-440	490	L	350-450	310-450
90° Standard	460	S	380-440	340-440	490	S	380-450	340-450	500	S	380-460	340-460
	490	S	390-450	350-450	500	S	390-460	350-460	510	S	390-470	350-470
	500	L	400-460	360-460	510	L	400-470	360-470	520	L	400-480	360-480
90° +50mm	500	S	400-460	360-460	510	S	400-470	360-470	520	S	400-480	360-480
	510	S	410-470	370-470	520	S	410-480	370-480	530	S	410-490	370-490
	520	L	420-480	380-480	530	L	420-490	380-490	540	L	420-500	380-500
75° tapered	430	S	370-450	NC	460	S	370-460	NC	500	S	370-470	NC
	460	S	380-460	NC	500	S	380-470	NC	510	S	380-480	NC
	500	L	390-470	NC	510	L	390-480	NC	520	L	390-490	NC
90° tapered	460	S	390-440	NC	490	S	390-450	NC	500	S	390-460	NC
	490	S	400-450	NC	500	S	400-460	NC	510	S	400-470	NC
	500	L	410-460	NC	510	L	410-470	NC	520	L	410-480	NC