

# ORDER FORM

## PIPER Comfort





Release date 2023/09



Company	Invoice address	Order number / reference	
Street		Customer reference	
E-Mail		Date	
City, Zip code		Offer	
Telephone		Order	
Company	Delivery address		
Street	City, Zip code		

Standard equipment:	Size / order codes
<ul style="list-style-type: none"> <li>• seat depth non-adjustable</li> <li>• extended backrest height non-adjustable</li> <li>• backrest recline adjustable</li> <li>• footrest height adjustable, detachable</li> <li>• seat width adjustable by FLEXI-pads set (optional eq.)</li> <li>• 5-point immo.harness with crotch belt and shoulder padding</li> <li>• 175/225 mm solid wheels, front with swivel lock</li> <li>• detachable wheels, forks</li> <li>• one-piece foldable handle</li> <li>• basket, transportation lock</li> </ul>	front swivel castor-wheels /FSCW
	PIC34PYYY €

Backrest height [mm]	680	with seat inlay / footrest raiser	680
Seat depth [mm]	290		240
Seat width [mm]	260/step/340		230/270
Seat to footrest [mm]	240-370		200/230
Backrest angle	90° - 125°		90° - 125°
Seat angle	25°		25°
Maximum load	50,0 kg		50,0 kg




colour frame / pattern colour	
frame	K2 - Anthracite
fabric patterns	G06  green
	G07  black
	G12  pink
	G13  blue

## PIPER Comfort Accessories

### Positioning and Fixation



#### Head

No.	Desc. / Codes	PIC34		
A1	Headrest adjustable RPRB026	UNI		€



#### Trunk

B1	SP - 5-point immobilisation belt RPRB009	UNI (SE)		€
B2	H-style harness RPRB007	S <sub>(1)</sub>	M	€



#### Pelvis / Abduction

C1	Abduction Block RPRB015	S <sub>(1)</sub>	M	€
C2	FLEXI-pads set (2H+2H) RPRB027	UNI		€
C3	Thigh fixation RPRB024	XS <sub>(1)</sub>	S	€
C4	Seat inlay- PIPER (seat size reduction) RPRB03001	UNI		€



#### Legs / Arms

No.	Desc. / Codes			
D1	Feet fixation - (pair) RPRB006	UNI		€
D2	Work board (soft) RPRB018	UNI		€
D3	Footplate booster (soft, 5x5cm height) RPRB031	UNI		€



D1



D2



D3



## Protection / transport

E1	Canopy Classic RPRB02901	UNI	€
E2	Raincover RPRB01701	M	€
E3	Footmuff - Summer RPRB028	UNI	€
E4	Footmuff - Winter RPRB025	M	€
E5	„TIE DOWN“ set - EU markets only RPRB03301	UNI	€
E6	„TIE DOWN“ set - USA, CDN markets only RPRB03302	UNI	€



E1




E2




E3



E4



E5



E6

### Sizes:

S - small

M - medium

UNI - universal size (one size only)

### Other notes:

(SE) Standard Equipment

(1) recommended with seat inlay or for smaller children