

ORDER FORM

CORZINO Classic

Release date 2023/08













patron[®]

creative in care

Company	Invoice address	Order number / reference	
Street		Customer reference	
E-Mail		Date	
City, Zip code		Offer	
Telephone		Order	
Company	Delivery address		
Street	City, Zip code		

Standard equipment:	Size / order codes			
	30	34	38	42
<ul style="list-style-type: none"> • 3D flexi footrest height / angle adjustable and foldable • backrest angle adjustable (2 positions) • 5-point immo. harness with crotch belt and shoulder padding • comfort side protection • backrest height extension (size 38, 42 only) • 175/265 mm solid wheels 	front swivel castor-wheels / FSCW			
	CNC30RP02	CNC34RP02	CNC38RP02	CNC42RP02
	€	€	€	€

Backrest height [mm]	550	650	650 (820)	660 (830)
Seat depth [mm]	280	330	380	380
Seat width [mm]	300	340	380	420
Seat to footrest [mm]	270-410	290-430	300-470	320-480
Backrest angle	91°, 93°	91°, 94°	90°, 93°	90°, 93°
Seat angle	22°, 21°	23°, 22°	22°, 21°	23°, 22°
Maximum load	45,0 kg	55,0 kg	75,0 kg	75,0 kg



colour frame / pattern colour							
frame	K2 - Anthracite						
fabric patterns	<table style="width: 100%; text-align: center;"> <tr> <td>G02  Tekno ruby red</td> <td>G03  Tekno dark blue</td> <td>G06  Tekno grass green</td> </tr> <tr> <td>G07  Tekno black</td> <td></td> <td></td> </tr> </table>	G02  Tekno ruby red	G03  Tekno dark blue	G06  Tekno grass green	G07  Tekno black		
	G02  Tekno ruby red	G03  Tekno dark blue	G06  Tekno grass green				
G07  Tekno black							

CORZINO Classic Accessories

Positioning and Fixation



Head

No.	Desc. / Codes	CNC 30	CNC 34	CNC 38	CNC 42	
A1	Headrest adjustable RPRB002	UNI	UNI	UNI	UNI	€
A2	Headwings BASIC RPRB013	M	M	L	L	€
A3	Backrest Height Extension			UNI(SE)	UNI(SE)	€




Trunk

B1	SP - 5-Point Immo. Belt CLASSIC RPRB010	S(SE)	M(SE)	M(SE)	M(SE)	€
B2	5-Point immo.belt COMFORT RPRB021	S	M	M	M	€
B3	Chest immo. belt RPRB004	S	M	M	L	€
B4	Immobilisation Waistcoat RPRB005	UNI	UNI	UNI	UNI	€
B5	H-Style Harness RPRB007	S	M	M	M	€
B6	Comfort side protection RPRB020	S (SE)	M (SE)	M (SE)	M (SE)	€
B7	Extended side protection RPRB032	S	M	L	L	€








Pelvis / abduction

No.	Desc. / Codes	CNC 30	CNC 34	CNC 38	CNC 42	
C1	Abduction Block RPRB015	S	M	M	M	€
C2	2-Point Safety Belt RPRB019	S	M	M	M	€
C3	Thigh fixation RPRB024	S	S	M	M	€





Legs / Arms

D1	Footrest 3D Flexi -2-piece angle adj. RPRB01130	UNI (SE)				€
D2	Footrest 3D Flexi, 2-piece angle adj. RPRB01134		UNI (SE)	UNI (SE)	UNI (SE)	€
D3	SP-Footrest 1-piece fold. OPTIONAL RPRB04030	UNI				€
D4	SP-Footrest 1-piece fold. OPTIONAL RPRB04034		UNI			€
D5	SP-Footrest 1-piece fold. OPTIONAL RPRB04038			UNI		€
D6	SP-Footrest 1-piece fold. OPTIONAL RPRB04042				UNI	€
D7	Feet Fixation - pair RPRB006	UNI	UNI	UNI	UNI	€
D8	Work board (soft) RPRB018	UNI	UNI	UNI	UNI	€
D9	Armrest foldable (set) RPRB02302	UNI	UNI			€
D10	Armrest foldable (set) RPRB02304			UNI	UNI	€



Protection / transport accessories

No.	Desc. / Codes	CNC 30	CNC 34	CNC 38	CNC 42	
E1	Canopy buggy CLASSIC RPRB00101	UNI	UNI	UNI	UNI	€
E2	Canopy buggy BASIC RPRB00102	UNI	UNI	UNI	UNI	€
E3	Rain cover - M RPRB01701	M	M			€
E4	Rain cover - L RPRB01702			L	L	€
E5	COMFORT Seat Cover RPRB003	M	M	L	L	€
E6	Basket - textile RPRB00801	UNI	UNI			€
E7	Basket - textile RPRB00802			UNI	UNI	€
E8	Footmuff - Summer RPRB012	M	M	L	L	€
E9	Footmuff - Winter RPRB025	S	M	M	L	€
E10	„TIE DOWN“ set - EU only RPRB01601	UNI(1)	UNI(1)	UNI(1)	UNI(1)	€
E11	„TIE DOWN“ set - USA, CDN only RPRB01602	UNI(2)	UNI(2)	UNI(2)	UNI(2)	€



Sizes:

XS - extra small
S - small
M - medium
L - large
XL - extra large
UNI - universal size (one size only)

Other notes:

(SE) Standard Equipment
(empty position) - accessory is not available
UNI(1) - offered for EU version only
UNI(2) - offered for USA and CDN version only